



Portable Hoist Connection System




















## Patented Connectus System as simple as 1,2,3...

- Using Connectus ballchain lower the strap to your desired height.
- 2. Connect portable hoist with swivel buckle to the Connectus strap and rewind strap to the ceiling.
- 3. Portable hoist is now ready for use. As no pole is in place, carers can move about freely and travel of hoist is not impeded.....



- 1300 789 420
- Sales@chshealthcare.com.au
- www.chshealthcare.com.au

## MELBOURNE

- Pallam, Victoria, 3083
- 1300 789 420

## SYDNEY

- 41 Edward Street Riverstone, NSW, 2765
  - 1300 789 420

## BRISBANE

- 7 Machinery Street Darra, Qld, 4076
- 1300 789 420









