

# L. DISPOSABLE SLING

With infection control becoming an area of greater focus and attention, Disposable Slings are designed to be used by one individual and disposed of when; a) the patient/resident no longer requires it or b) the sling is soiled.

The Disposable Slings are based on the design, safety, comfort and application of the Deluxe Universal Sling pattern. All disposable slings are rated to 272kg. Label indicates when the sling is in need of disposal. Typical safe use ranges from 1-3 months (dependant upon frequency of use).

**For use with standard yoke and four point yoke.**

**Colour: White**

|                   | S      | M      | L      | XL     |
|-------------------|--------|--------|--------|--------|
| A (mm)            | 940    | 1065   | 1195   | 1320   |
| B (mm)            | 355    | 355    | 355    | 355    |
| C (mm)            | 585    | 650    | 710    | 775    |
| D (mm)            | 595    | 675    | 815    | 950    |
| E (mm)            | 1535   | 1675   | 1880   | 2085   |
| F (mm)            | 470    | 470    | 470    | 470    |
| G (mm)            | 200    | 255    | 305    | 355    |
| H (mm)            | 255    | 255    | 255    | 255    |
| PRODUCT CODE-10pk | 547050 | 547060 | 547070 | 547080 |
| SWL (Kg)          | 272    | 272    | 272    | 272    |

