

Creating the optimal work environment with TurnAid T4

Main uses of TurnAid system

TurnAid enhances and streamlines the care, relocation and accommodation processes for individuals with severe disabilities in care facilities, institutions, and private residences.

How does the TurnAid system work

TurnAid is an assistive advice installed on care beds. The system integrates two established core functions: a siderail and an electric turning mechanism. TurnAid is easily installed on most care bed models, thanks to its universal brackets and adjustable length and width. The system complies with the rigorous standards of the associated medical directive.

TurnAid accommodates individuals weighing up to 350kg, with a focus on ensuring the safety and comfort of both the user and caregiver throughout the task.



TurnAid's key features

og mulighed for ekstra tilbehør

TurnAid comes equipped with

- 1 set of siderails featuring an integrated electric turning system (both right and left and sides)
- 1 installation cover with electric parts:
 1 control box, 2 hand controls and assorted engine cables
- 1 foot control
- 2 slide sheets
- · 2 pull sheets



Foot control

Item no: 3010 HMI: 125694



Hand control

Item no: 3030

TurnAid Accessories - optional to buy

- Foot control*
- Hand control*
- · Siderail guards
- Slide sheets available in six different sizes
- Pull sheets available in three different sizes
- Transport and assembly trolley
- Spacers
- Mounting bolts

*for hand and foot controls the speed of the motor rollers can be adjusted to 100%, 75% and 50% using the Speed Control function.



Siderail guards

Item no: 5050 HMI: 89929



Slide sheet

Item no: 2201 size: 95 x 205 cm HMI: 89879



Pull sheet

Item no: 2552 size: 180x350 cm HMI: 89881



Transport and assembly trolley

Item no: 6003 HMI: 131243



Spacers

Item no: 4093



Mounting bolts 60 mm

Item no: 4096

We aim to provide the best workplace solutions



With both the individual and the caregiver in mind, we deliver the highest standards of care – for this reason we are also ISO certified.

We continuously strive to promote compassion and comfort in both care and handling. With a steadfast commitment to customer satisfaction, we aim to become your partner of choice.

Danrehab A/S
Vejlevej 15 | Ølholm | DK-7160 Tørring
Tel. +45 86 44 73 55

info@danrehab.dk danrehab.dk

